



QLABTICS

SmileFuge Micro Centrifuge



Register in [www.qlabtics.com](http://www qlabtics.com)
then scan me



Qlabtics LLC
CHICAGO RIDGE, IL 60415

SmileFuge, a Micro Centrifuge, comes with a free QR code for equipment maintenance and calibration in labs. Qlabtics provides a range of cloud-based solutions to manufacturers, service providers, and lab directors to improve equipment management, reduce downtime, and increase customer satisfaction. Their solutions include equipment tracking and maintenance, calibration management, remote equipment monitoring, engineer management, customer follow-up, asset tracking, live device usage tracking, cloud-based reporting, and document management. Visit their website at www.qlabtics.com.



SmileFuge Micro Centrifuge



The Smilefuge microcentrifuge is the perfect choice for quick spin-down applications and microfiltration. With easy-to-use operation, it begins rotating as soon as the lid is closed. It is compatible with both 2mLx8 and PCR8x4 rotors and runs quietly with a noise level of $\leq 45\text{Db}$. Changing the rotor is a breeze, requiring no tools, and acceleration and braking take only seconds, durable rotor with a clamp locking design that ensures safety.

Features

- The Smilefuge microcentrifuge is ideal for microfiltration and quick spin-down applications.
- It is easy to operate, starting to rotate when the lid is closed.
- The centrifuge is compatible with 2mLx8 and PCR8x4 rotors.
- The rotor can be quickly changed without any tools.
- The centrifuge accelerates and brakes in seconds.
- It is available in two colors, green and blue.
- The rotor has a clamp locking design for safety and durability.



SmileFuge	QSF-007
Max. Speed	7000rpm
Rotor Capacity	0.2/0.5/1.5/2.0mLx8; 32x0.2mL PCR strips or 4x0.2mL PCR 8 strips
Run Time	Continuous operation
Power	AC100V-240V/50Hz/60Hz 20W
Noise Level	$\leq 45\text{dB}$
Dimension [DxWxH]/ weight	160mmx170mmx122mm/ 0.5Kg

Website: us.qlabtics.com
Email: info@qlabtics.co

QLABTICS LLC
10500 PARKSIDE AVE
APT 9 CHICAGO RIDGE,
IL 60415 USA

Register in www.qlabtics.com
then scan me



Qlabtics LLC
CHICAGO RIDGE, IL 60415