

## Accessories for VortexBlind Vortex Mixers

### VT1.3.1 Test tube adapter Article number: 18900020

48 × Ø 6mm (1mL test tube)



### VT1.3.2 Test tube adapter Article number: 18900021

15 × Ø 10mm (1.5 - 2mL test tube)



### VT1.3.3 Test tube adapter Article number: 18900022

16 × Ø 12mm (5mL test tube)



### VT1.3.4 Test tube adapter Article number: 18900023

8 × Ø 16mm (10mL test tube)



### VT1.3.5 Test tube adapter Article number: 18900024

8 × Ø 20mm (10 - 15mL test tube)



### VT1.3.6 Flat gasket Article number: 18900043

Ø 99mm, suitable for multi tube jog



### VT1.3.7 Vacuum suction cup feet Article number: 18900158

Material: Rubber



### VT1.3.8 Test tube adapter Article number: 18900395

24 well, Ø 10mm (1.5-2mL test tube)



### VT1.3.9 PCR Octuple Tube Adapter Article number: 18204185

56 × Ø 5mm (0.1 - 0.5mL test tube)



### VT1.3.10 test tube adapter Article number: 18900855

4-well, Ø 28mm (50mL test tube)



### VT1.3.11 Multi specification test tube adapter Article number: 18900862

12×Ø 9.8mm+12×Ø 6.6mm+12×Ø 5mm (1.5/2.0/0.5/0.2mL test tube)



### VT1.2 adapter wrench Article number: 18900044

Used for loading and unloading vibration head adapters



**VT1.3 Universal Tray Article number: 18900505**

Ø 100mm



**PS1.0 single hole plate tray Article number: 18204242**

Work board (PCR board with full skirt edge, cell culture board, enzyme-linked immunosorbent assay board, deep well plate, etc., half skirt edge/no skirt edge PCR board needs to be used with well plate adapter)



**VT1.3.12 large platform test tube adapter Article number: 18900875**

44×Ø 5mm+55×Ø 6.6mm+30×Ø 9.8mm(1.5/2.0/0.5/0.2mL)



**VT1.4 Ampoule/Test Tube Accessories Article number: 18204566**

4 ampoule bottles/test tubes (suitable for 0.2mL-15mL test tubes)



**VT1.5 Conical Bottle Attachment Article number: 18204561**

Suitable for 100mL-250mL conical flask, recommended rotation speed: 0-1500rpm

