



QLABTICS

MultiSpin Mini Centrifuge 3-10K



Register in www.qlabtics.com
then scan me



Qlabtics LLC
CHICAGO RIDGE, IL 60415

MultiSpin Mini Centrifuge 3-10K, comes with a free QR code for equipment maintenance and calibration in labs. Qlabtics provides a range of cloud-based solutions to manufacturers, service providers, and lab directors to improve equipment management, reduce downtime, and increase customer satisfaction. Their solutions include equipment tracking and maintenance, calibration management, remote equipment monitoring, engineer management, customer follow-up, asset tracking, live device usage tracking, cloud-based reporting, and document management. Visit their website at www.qlabtics.com.

MultiSpin Mini Centrifuge 3-10K

Introduction:

the MultiSpin Mini Centrifuge, a sleek and innovative centrifuge designed to meet diverse laboratory needs. With its unique appearance and user-friendly features, this centrifuge offers exceptional performance and convenience.

Features:

- Speed adjustable from 3000-10000rpm.
- Equipped with two types of rotors and multiple tube carriers.
- Compatible with 2.0ml, 0.5ml, 0.2ml tubes, and 0.2ml x 8 PCR tube strips.
- 32 - holes Rotor.
- Built-in timing function for precise experimental control.
- Large radius rotor with 8 tube holes for increased centrifugal force.
- Press-button cover opens automatically to 95° for easy one-hand operation.
- Clip-on fixed rotor for convenient and reliable rotor exchanges.
- LED displays for speed and timing.
- Silent and stable operation with a maintenance-free motor,
- Versatile design to meet various experimental demands.



Website: us.qlabtics.com
Email: info@qlabtics.co

QLABTICS LLC
10500 PARKSIDE AVE
APT 9 CHICAGO RIDGE,
IL 60415 USA

Register in www.qlabtics.com
then scan me



Qlabtics LLC
CHICAGO RIDGE, IL 60415

AquaTemp Thermostatic Water Bath



Model	QMS-010
Max. Speed	3000-10000rpm
Max. RCF	500-5400xg
Rotor	(16) 2.0/1.5/0.5/0.2ml centrifuge tubes (With 0.5ml and 0.2ml tube adapter) (16) 0.2ml PCR Strips(when speed≤6000rpm) 32 - hole Rotor.
Noise	≤55dB
Power Supply	AC110/220V±10%, 50/60Hz
External Size(W*D*H)mm	180*160*121
Package Size(W*D*H)mm	280*220*200
Net/Gross Weight(kg)	1.4/3

***2 years warranty.**

**** free of charge, unique software of Qlabtics for managing device maintenance and troubleshooting.**

Website: us.qlabtics.com
Email: info@qlabtics.co

QLABTICS LLC
10500 PARKSIDE AVE
APT 9 CHICAGO RIDGE,
IL 60415 USA

Register in www.qlabtics.com
then scan me



Qlabtics LLC
CHICAGO RIDGE, IL 60415