



QLABTICS

TubeFlexor Tube Roller



Register in [www.qlabtics.com](http://www qlabtics.com)
then scan me



Qlabtics LLC
CHICAGO RIDGE, IL 60415

TubeFlexor Tube Roller, comes with a free QR code for equipment maintenance and calibration in labs. Qlabtics provides a range of cloud-based solutions to manufacturers, service providers, and lab directors to improve equipment management, reduce downtime, and increase customer satisfaction. Their solutions include equipment tracking and maintenance, calibration management, remote equipment monitoring, engineer management, customer follow-up, asset tracking, live device usage tracking, cloud-based reporting, and document management. Visit their website at www.qlabtics.com.

TubeFlexor Tube Roller

TubeFlexor by Qlabtics is designed for uniform mixing of viscous liquids, whole blood components, and blood samples contained in anticoagulation tubes or bottles. This versatile tube roller can mix samples in any type of liquid container. Its microprocessor control and digital timer ensure precise and reliable mixing every time.

Features:

- + Maximum load of up to 4Kg.
- + Gentle mixing in rolling motion with a capacity of 6 tubes.
- + Speed adjustable from 0-70rpm.
- + LED overload indicator



Features:

- + **Maximum load of up to 4Kg.**
- + **Gentle mixing in rolling motion with a capacity of 6 tubes.**
- + **Speed adjustable from 0-70rpm.**
- + **LED overload indicator**

*2 years warranty.

** free of charge, unique software of Qlabtics for managing device maintenance and troubleshooting.

Website: us.qlabtics.com
Email: info@qlabtics.co

QLABTICS LLC
10500 PARKSIDE AVE
APT 9 CHICAGO RIDGE,
IL 60415 USA

Register in www.qlabtics.com
then scan me



Qlabtics LLC
CHICAGO RIDGE, IL 60415

TubeFlexor Tube Roller



Specifications	TF-001
Amplitude	24mm
Motor type	DC motor
Rocking motion	Rocking and rolling
Amplitude	24mm
Max. load capacity	4kg
Number of rollers	6
Roller size[length]	280mm
Speed range	0 -70rpm
Speed display	Scale
Operation mode	Continuous
Voltage	100-240V,50/60Hz
Power	25W
Dimension [D×W×H]	450×260×120mm
Weight	4.5kg
Permissible ambient temperature and humidity	5-40°C,80%RH
Protection class	IP21

Website: us.qlabtics.com
Email: info@qlabtics.co

QLABTICS LLC
10500 PARKSIDE AVE
APT 9 CHICAGO RIDGE,
IL 60415 USA

Register in www.qlabtics.com
then scan me



Qlabtics LLC
CHICAGO RIDGE, IL 60415