



VERTICAL PRESSURE STEAM STERILIZER

Characteristics:

- 1.Fully stainless steel structure
- 2.Indicator light indicates working state
- 3.Over temperature&over pressure auto-protection
- 4.Double scale indication pressure gauge
- 5.Safe protection of water lacking
- 6.Automatically shut off with beep reminding after sterilization
- 7.Easy to operate,safe and reliable
- 8.With two stainless steel sterilizing baskets

Technical Data:

Model technical data	LS-35LJ	LS-50LJ	LS-75LJ	LS-100LJ
Chamber volume	35L(ϕ 318×450) mm	50L(ϕ 340×550) mm	75L(ϕ 400×600) mm	100L(ϕ 440×650) mm
working pressure	0.22MPa			
Working temperature	134℃			
Max working pressure	0.23 Mpa			
Heat average	$\pm 1^{\circ}\text{C}$			
Timer scope	0~60min			
Temperature scope	105~134℃			
Power	2.5KW/AC220V.50Hz	3KW /AC220V.50Hz	4.5KW /AC220V.50Hz	
Overall dimension	470×450×840(mm)	520×500×970(mm)	550×530×990(mm)	580×560×1060(mm)
Transport dimension	540×540×970(mm)	580×580×1090(mm)	610×610×1110(mm)	640×640×1170(mm)
G.W/N.W	61kg/ 46Kg	73kg/56Kg	87kg/67Kg	95kg/73kg